




SIMPLE FLYWHEEL

HTM9



Year 1
study

Features

- Wall mounted apparatus.
- Can be mounted to HST1, HST100, HAC10, HAC14, HVT12f frames
- Flywheel 'locked' before weights applied
- All hangers and weights supplied
- Timing system supplied

Description

A steel disc Ø250mm diameter and 30mm thick is mounted on a shaft running in precision bearings housed in a substantial wall bracket. A mark on the flywheel and a pointer on the bracket enables the revolutions to be counted and timed with the stop watch supplied. A cord, Load hanger and set of weights are provided to start the flywheel rotating. From the timings of a known number of revolutions, students can then verify the experimental capabilities of the apparatus.

Related Laws/Applications

- Newton's Second Law of Motion
- Rotational Mass
- Rotational Acceleration
- Moment of Inertia
- Angular Acceleration
- Torque
- Couple
- Rotating Machinery

Learning capabilities

- To verify the second law of motion applied to a flywheel, i.e. the relationship between torque and angular acceleration
- To compare experimental and calculated moments of inertia of a disc
- To study the energy transformations and to demonstrate that a flywheel can be used to store energy

Technical Specification

- Coated Steel Flywheel
- Flywheel: Ø250 x 30mm thick
- Weights set: 10 x 0.1N, 10 x 0.2N, 1 x 1N, 2 x 2N, 1 x 5N, 1 x 10N
- 1 x Load hanger

Essential Ancillaries

- HAC14
- HAC10
- HST1
- HST100
- HVT12f

What's in the Box?

- 1 x HTM9 assembled
- 1 x Load hanger
- 10 x 0.1N, 10 x 0.2N, 1 x 1N weights
- 2 x 2N, 1 x 5N, 1 x 10N weights
- 1 x Stopwatch
- 1 x Tape measure
- Spare Cord
- Instruction manual
- Packing list
- Test sheet

Weights & Dimensions

- Weight: 18 kg
- Length: 190mm
- Width: 250mm
- Height: 300mm

Essential Services

- Recommended Suspension height not to exceed 1.0metres

Operational Conditions

- Storage temperature: -10°C to +70°C
- Operating temperature range: +10°C to +50°C
- Operating relative humidity range: 0 to 95%, non condensing

Ordering information

To order this product, please call PA Hilton quoting the following code:
HTM9

All brand and/or product names are trademarks of their respective owners. Specifications and external appearance are subject to change without notice. The colour of the actual product may vary from the colour shown in the brochure.

Copyright © 2018 P.A. Hilton Limited. All rights reserved. This technical leaflet, its contents and/or layout may not be modified and/or adapted, copied in part or in whole and/or incorporated into other works without the prior written permission of P. A. Hilton Limited. Hi-Tech Education is a registered trade mark of P. A. Hilton Limited.

COUNTRY OF ORIGIN - UK WARRANTY PERIOD - 5 YEARS