

## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled

Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils


Aluminium Anodized

Temperature ranges 0 - +130 °C

*(Other temperature on request)*

| Part Number | Description                                         | Sample Thickness | Rating | Fixing   | Weight |
|-------------|-----------------------------------------------------|------------------|--------|----------|--------|
| 01/4184     | Eccentric roller grip with crosscut faces 50mm wide | 0 - 7mm          | 1kN    | 5/8 inch | 509g   |

| Dimensions |                         |
|------------|-------------------------|
| Width      | 50mm                    |
| Height     | 80mm                    |
| Depth      | 76mm (including handle) |



## Eccentric Roller Grips



A range of self tightening eccentric roller grips, the grip tightens as the sample is pulled  
Ideal for all types of flat samples, rubber, plastic film, paper and cardboard and foils

Aluminium Anodized

Temperature ranges 0 - +70 °C (restricted due to rubber faces)

| Part Number | Description                                                                  | Sample Thickness | Rating | Fixing   | Weight |
|-------------|------------------------------------------------------------------------------|------------------|--------|----------|--------|
| 01/4185     | Eccentric roller grip with plain roller and rubber lined grip body 50mm wide | 0 – 7mm          | 1Kn    | 5/8 inch | 509g   |

| Dimensions |                         |
|------------|-------------------------|
| Width      | 50mm                    |
| Height     | 80mm                    |
| Depth      | 76mm (including handle) |





## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled  
Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils

Aluminium Anodized  
Temperature ranges 0 - +130 °C  
(Other temperature on request)

| Part Number | Description                                          | Sample Thickness | Rating | Fixing   | Weight |
|-------------|------------------------------------------------------|------------------|--------|----------|--------|
| 01/4186     | Eccentric roller grip with crosscut faces 100mm wide | 0 – 7mm          | 5kN    | 5/8 inch | 500g   |

| Dimensions |                         |
|------------|-------------------------|
| Width      | 100mm                   |
| Height     | 96mm                    |
| Depth      | 89mm (including handle) |



## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled

Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils


Aluminium Anodized

Temperature ranges 0 - +130 °C

(Other temperature on request)

| Part Number | Description                                         | Sample Thickness | Rating | Fixing   | Weight |
|-------------|-----------------------------------------------------|------------------|--------|----------|--------|
| 01/4187     | Eccentric roller grip with crosscut faces 50mm wide | 0 – 7mm          | 5kN    | 5/8 inch | 667g   |

| Dimensions |                         |
|------------|-------------------------|
| Width      | 50mm                    |
| Height     | 96mm                    |
| Depth      | 91mm (including handle) |



## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled

Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils


Aluminium Anodized

Temperature ranges 0 - +130 °C

(Other temperature on request)

| Part Number | Description                                         | Sample Thickness | Rating | Fixing   | Weight |
|-------------|-----------------------------------------------------|------------------|--------|----------|--------|
| 01/4188     | Eccentric roller grip with crosscut faces 60mm wide | 0 – 8mm          | 10kN   | 5/8 inch | 900g   |

| Dimensions |                          |
|------------|--------------------------|
| Width      | 90mm                     |
| Height     | 157mm                    |
| Depth      | 126mm (including handle) |



## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled

Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils

The 01/4189 is fitted with a sample alignment device to ensure samples are always fitted correctly by the operator


Aluminium Anodized

Temperature ranges 0 - +130 °C

*(Other temperature on request)*

| Part Number | Description                                                                     | Sample Thickness | Rating | Fixing   | Weight |
|-------------|---------------------------------------------------------------------------------|------------------|--------|----------|--------|
| 01/4189     | Eccentric roller grip with crosscut faces 60mm wide and sample centering device | 0 – 8mm          | 10kN   | 5/8 inch | 900g   |

| Dimensions |                          |
|------------|--------------------------|
| Width      | 90mm                     |
| Height     | 157mm                    |
| Depth      | 126mm (including handle) |



## Eccentric Roller Grips

A self tightening eccentric roller grips, the grip tightens as the sample is pulled  
Ideal for all types of flat samples, rubber, plastic film, plastic dumbbells, paper and cardboard, sheet metal and foils  
The 01/4190 has dual rollers for symmetrical gripping of the sample

Aluminium Anodized  
Temperature ranges 0 - +130 °C  
(Other temperature on request)

| Part Number | Description                            | Sample Thickness | Rating | Fixing   | Weight |
|-------------|----------------------------------------|------------------|--------|----------|--------|
| 01/4190     | Eccentric dual roller grip 30mm x 50mm | 0 – 7mm          | 5kN    | 5/8 inch | 1.1Kg  |

