

# Thread Head Grips

Thread head grips consist of a two part fixture, the main body is chosen to fit the eye end of the machine and Loadcell and the insert is chosen to fit the diameter of the thread on the sample. Samples are machined dumbbells with metric threads

Inserts are also available for bolts and screws

Steel  
Temperature ranges -70 - +230 °C  
*(Other temperature on request)*

Part Number	Description	Rating	Fixing
01/4717	15.9mm (5/8 inch) eye end adaptor	20kN	5/8 inch
01/4718	Inserts, M3, M4, M5, M6 and M8	20kN	5/8 inch
01/4719	31.8mm (1 ¼ inch) eye end adaptor	100kN	1 ¼ inch
01/4720	Inserts in M3 to M16 in 1mm intervals	100kN	1 ¼ inch



